

Hemodynamic Monitoring: The Trauma Patient

Clinical Profile:

Trauma (traw'mă) [Gr. trauma, wound]. 1. A physical injury or wound caused by external force or violence. 2. An emotional or psychological shock that may produce disordered feelings or behavior.¹

Traumatized patients with bodily injury or illness require expedient, thorough and advanced critical healthcare. Essential, therefore, is the vital accuracy in hemodynamic measurement of preload to optimize cardiac output and oxygen delivery.² Such assessment has historically been performed using filling pressures [central venous pressure (**CVP**) and pulmonary artery wedge pressure (**PAWP**)], due to an assumed correlation with end diastolic volume.³ A new generation of volumetric measurements, right ventricular ejection fraction (**RVEF**) and right ventricular end diastolic volume (**RVEDV**), however, offers another method of assessment, which is evaluated herein.

Correlations:

• Correlation #1: **RVEDVI** and **SI / CI**

Following volume administration over 30 minutes, changes in **SI** have been shown not to correlate with **RAP** or **PAWP** but did correlate with **RVEDVI**. In patients with **RVEDVI** greater than 140 ml/m², the fluid challenge increased **RAP** and **PAWP**, and decreased **LVSWI**.⁴ Discrepancies between **PAWP** and **RVEDVI** occurred in more than 50% of patients studied and demonstrated that **RVEDVI** more accurately predicted preload recruitable increases in **CO**.² In a study conducted by Karen Safcsak, R.N., et al., it was surmised that the **RVEDVI** (rather than pulmonary artery occlusion pressure [**PAOP**]) better correlates with **CI** in assessing the ventricular preload of surgical, trauma, post shock volume-resuscitated and respiratory failure patients, making **RVEDVI** "the new gold standard in the assessment of ventricular preload".⁵ Studies have confirmed the superiority of **RVEDVI** over **RAP** as an indicator of **RV** preload.⁶

• Correlation #2: **RVEF** and trauma survival

In trauma patients, survival has been shown to depend upon right ventricular function. In a study of 17 trauma patients, patients demonstrated a decrease in **RVEF**. In survivors, **RVEF** improved over 8-12 hours; in non-survivors, **RV** function continued to deteriorate.⁷

• Correlation #3: **RVEDVI** and trauma survival, intestinal perfusion and resuscitation

The results of a resuscitation study by Dr. Michael Cheatham, et al., demonstrated that optimal oxygen delivery and tissue perfusion, regardless of ventricular function, is a more individualized (per patient) method for resuscitation, than utilizing **RVEDVI**, **CVP**, or **PAOP** measurements alone for volume resuscitation.⁹ Trauma patients resuscitated to an **RVEDVI** of 120 ml/m² have significantly better visceral perfusion than those patients resuscitated to 90 to 100 ml/m² with inotropic augmentation of blood pressure as needed.⁸

A study by Dr. Chang, et al., showed that **RVEDVI** of >100ml/m² = ↑ intestinal perfusion, ↓ mortality, ↓ MOF.¹⁰

Outcome:

It has been recognized that a new catheter, combining **SvO₂** with **RVEF**, **RVEDV** measurements and continuous thermodilution cardiac output (**CTCO**), will yield not only information regarding oxygen consumption and delivery, but will also permit a better understanding of hemodynamics. Loren D. Nelson, M.D. concluded in his report that the combination of intravascular pressure measurements and continuous cardiac output (**CCO**) also allows calculation of left ventricular (**LV**) and right ventricular (**RV**) stroke work that, "under conditions of steady state preload and afterload, reflect the contractile function of the heart as a determinant of cardiac performance." On a near continuous basis, the critical care team will thus be provided with complete hemodynamics and oxygen transport values.³



Edwards
LIFESCIENCES

The quotes and information referenced in this brochure were selected from independent third party publications and are not intended to suggest that such third parties have reviewed and/or endorsed Edwards' products.

References:

1. Taber, CW (1977): Taber's Cyclopedic Medical Dictionary, (13th edition): Philadelphia: FA Davis Company
2. Diebel LN, et al., "End-Diastolic Volume: A Better Indicator of Preload in the Critically Ill", Arch Surg, 1992;127:817-822.
3. Nelson LD., "The New Pulmonary Arterial Catheters: Right Ventricular Ejection Fraction and Continuous Cardiac Output", Critical Care Clinics, 1996;12:795-818.
4. Reuse C, et al., "Measurements of Right Ventricular Volumes During Fluid Challenge." Chest, 1990;98:1450-1454.
5. Safcsak K, et al., "Right Heart Volumetric Monitoring: Measuring Preload in the Critically Injured Patient", AACN Clinical Issues, 1999;10:22-31.
6. Edwards JD., "Time for a Reappraisal of the Use of the Right Ventricular Ejection Fraction Thermodilution Catheters?", Critical Care Medicine, 1998;26:1769-1770.
7. Eddy AC, et al., "Right Ventricular Dysfunction in Multiple Trauma Victims", Am J Surg, 1988; 155:712-715.
8. Miller PR, et al., "Randomized Prospective Comparison of Increased Preload Versus Inotropes in the Resuscitation of Trauma Patients: Effects on Cardiopulmonary Function and Visceral Perfusion", J. Trauma, 1998;44:107-113.
9. Cheatham ML, et al., "Right Ventricular End-diastolic Volume Index is Superior to Cardiac Filling Pressures in Determining Preload Status", Surg Forum, 1994;45:79-81.
10. Chang MC, et al., "Cardiac Preload, Splanchnic Perfusion, and their Relationship During Resuscitation in Trauma Patients", J Trauma, 1997; 42: 577-584.



**Certified according
to European MDD**

See instructions for use for full prescribing information.

Edwards Lifesciences, Edwards, and the stylized E logo are trademarks of Edwards Lifesciences Corporation.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

© Copyright 2003 Edwards Lifesciences LLC
All rights reserved. 3163A-12/01-CC



Edwards Lifesciences

Edwards Lifesciences LLC · One Edwards Way · Irvine, CA 92614 USA · 949.250.2500 · 800.424.3278 · www.edwards.com
Edwards Lifesciences (Canada) Inc. · 1290 Central Pkwy West, Suite 300 · Mississauga, Ontario · Canada L5C 4R3 · Phone 905.566.4220 · 800.268.3993
Edwards Lifesciences S.A. · Ch. du Glapin 6 · 1162 Saint-Prex · Switzerland · Phone 41.21.823.4300
Edwards Lifesciences Japan · 2-8 Rokubancho · Chiyoda-ku, Tokyo 102-0085 · Japan · Phone 81.3.5213.5700