


PediaSat Oximetry Catheter

The PediaSat oximetry catheter is the first oximetry catheter designed for pediatrics that can provide early indication of critical changes in a patient's clinical condition via continuous, real-time monitoring of oxygen delivery and consumption.⁴⁻⁵

It helps guide therapy and enables early intervention in your most complex pediatric patients. The PediaSat oximetry catheter is designed for use with the Vigileo monitor.



PediaSat oximetry catheter* specifications


Model Number	Lumens	Length cm	Size F (mm)
XT245HS	2	5	4.5 (1.50)
XT248HS	2	8	4.5 (1.50)
XT358HS	3	8	5.5 (1.83)
XT3515HS	3	15	5.5 (1.83)
The following kits are available in the U.S. only.			
XT245HK	2	5	4.5 (1.50)
XT245K	2	5	4.5 (1.50)
XT248HK	2	8	4.5 (1.50)
XT248K	2	8	4.5 (1.50)
XT358HK	3	8	5.5 (1.83)
XT358K	3	8	5.5 (1.83)
XT3515HK	3	15	5.5 (1.83)
XT3515K	3	15	5.5 (1.83)

* PediaSat catheters are designed for use with Edwards oximetry monitors and optical cables. ** All model numbers with an "H" contain AMC Thromboshield, an antibacterial heparin coating which decreases visible microbe count on the surface of product during handling and placement. Many catheters are available with or without heparin coating.

Continuous ScvO₂ Monitoring with the PediaSat Oximetry Catheter

Continuously monitoring central venous oxygen saturation (ScvO₂) enables the early detection and management of occult tissue hypoxia.^{1,2}

Components of ScvO₂



- Reinhart, K, et al. Continuous central venous and pulmonary artery oxygen saturation monitoring in the critically ill. *Intensive Care Med* 2004;30(8):1572-8.
- Zaja, J., Venous oximetry. *SIGNA VITAE* 2007; 2(1):6-10
- Loren, D. Continuous Venous Oximetry in Surgical Patients. *Ann Surg* 1986;203/3:329-333.
- Krahn, G, et al. Continuous central venous saturations during pericardial tamponade case report. *Pediatr Crit Care Med* 2007, Vol. 8, No. 3, p18.2.153.
- Krahn, G, et al. Early clinical evaluation of the Edwards PediaSat™ oximetry catheter in pediatric patients. *Pediatr Crit Care Med* 2007, Vol. 8, No. 3, p18.2.152.

Rx only. See instructions for use for full prescribing information. Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity. Edwards is a trademark of Edwards Lifesciences Corporation. Edwards Lifesciences, the stylized E logo, PediaSat and Vigileo are trademarks of Edwards Lifesciences Corporation and are registered in the United States Patent and Trademark office. ©2008 Edwards Lifesciences LLC. All rights reserved. AR03544



Edwards Lifesciences LLC · Irvine, CA 92614 USA · 949.250.2500 · www.edwards.com
 Edwards Lifesciences (Canada) Inc. · Mississauga, Ontario · Canada L5C 4R3 · 905.566.4220
 Edwards Lifesciences Europe · 1162 Saint-Prex · Switzerland · 41.21.823.4300
 Edwards Lifesciences Japan · Chiyoda-ku, Tokyo 102-0085 · Japan · 81.3.5213.5700

